ZIU04965

INTERPUPILLOMETRO DIGITALE

DIGITAL PD METER / PD-MESSGERÄT / APPAREIL DE MESURE DU DP

Ref. **04965**

MANUALE D'USO E MANUTENZIONE - GARANZIA
OPERATION AND MAINTENANCE MANUAL - WARRANTY
GEBRAUCHS - UND WARTUNGSHANDBUCH
MANUEL D'INSTRUCTIONS ET D'ENTRETIEN - GARANTIE

OPTICS & SERVICE
centrostyle®

CENTROSTYLE S.p.A. - Via G.D.Martinengo,7 - 21040 Veduggio Olona - Varese - Italy

Servizio Clienti Italia: Tel. 0332.270270 - Fax 0332.270222 E-mail: infocli@centrostyle.it

Customer Service Export: Tel. +39.0332.270400 - Fax +39.0332.270405 E-mail: export.dtp@centrostyle.it

Kundenservice Italien: Tel. +39.0332.270400 - Fax +39.0332.270405 E-mail: export.dtp@centrostyle.it

Service Clients Italie: Tel. +39.0332.270400 - Fax +39.0332.270405 E-mail: export.dtp@centrostyle.it

www.centrostyle.com

SERVICE CENTRE

OPTICS & SERVICE
centrostyle®

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Servizio clienti: Tel. 0332 270270 (orario: 8.30/12.30 - 14.00/18.00)

Page 1	TRANSPORT AND STORAGE.
Page 1	STANDARD ACCESSORIES.
Page 2	DESCRIPTION OF THE INSTRUMENT.
Page 2	SETTING UP, START-UP AND OPERATION.
Page 6	SAFETY PRECAUTION.
Page 6	MAINTENANCE AND CLEANING.
Page 7	PUTTING OUT OF SERVICE.
Page 7	TECHNICAL SPECIFICATION.
Page 8	WARRANTY AND SERVICE.

TRANSPORT AND STORAGE ---

Always ship and store the instrument and accessories in their original packaging.

- Do not expose the instrument to temperatures below 0° or above 60° C, to relative humidity above 85% or to a very dusty surrounding.
- Do not stock more than 15 packages of weights in excess of 13 kg on top of an instrument's packaging. Do not subject to hard knocks or mechanical strain.

STANDARD ACCESSORIES ---

1 digital PD-meter ref. 04965
1 operation manual ref. 81015
4 batteries 1,5 V AA (LR6)

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DESCRIPTION OF THE INSTRUMENT

This instrument is a digital PD-meter, both for binocular and monocular measurements. It is an ancillary instrument in the determination of the optical centre of ophthalmic lenses for correction glasses.

TARGET USERS: professional, only for use by adequately trained staff.

APPLICATION FIELD: optical shops, optical laboratories. Instrument for external measurement, not a medical diagnostic instrument.

SETTING UP, START-UP AND OPERATION


Open the packaging and make sure it is undamaged, check that it contains all parts listed on page 1 and are in perfect condition.

If damage has incurred during transport, contact Centro Style or its distributor.

- Open both battery compartments and insert the 4 batteries included in the correct direction.


Description of all information and icons appearing on the front panel in the various modalities:


1. Activate the instrument by pressing reset key ②: The default PD (63 mm) appears ①, together with the relative right and left single PDs. The inside reference lamp lights up.

2. Set the instrument's nose bearer ⑫ lightly on the testee's nose so that their eyes correspond to the patient's window ⑬. The front bearer ⑨ must touch the testee's front.
3. The testee must continuously look at the green reference light inside the instrument.

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


4. The tester, through the measuring window ⑥ must align the testee's pupil with the reference lines in the patient's window ⑫

- The reference lines can be moved by pressing the PD measuring keys ③ and ④. Pressing the keys' inside, the PD diminishes, pressing the outside, the PD increases.
- Note: when the reference line reaches the lower or higher limit, the end darkens to advise the tester that the instrument is out of range.

5. Remove the instrument from the testee and note the total PD as well as the right and left PD.
6. Adapt the PD in relation to the convergence distance required for the spectacles in question. For this operation it is not necessary to measure again the PD on the testee as long as you do not press keys ③ or ④:

- To obtain the new PD press as many times the TBIS key ⑤ until the display shows the required distance. The convergence range goes from 30 cm to infinity.
- Read the new PD shown on the display.

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
- **Reset:** pressing the key ② the instrument can be rest in any moment to the default PD (63 mm), binocular modus and infinite convergence
- **Memory function:** pressing the key MEM ⑧ at the end of measurement, the data on the display is locked even if by error you press keys ③ or ④. To unlock press key MEM ⑧ or reset ②.
- **Monocular or binocular measurement:** press key BIN ⑦ to change to:
 - monocular left PD (the right ocular dims)
 - monocular right PD (the left ocular dims)
 - total PD value

- During measurement it is recommended to lower the baffles ⑪ to separat the tester and the testee's face

After 1 minute of non utilization, the instrument shuts down automatically.

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SAFETY PRECAUTIONS

- Replace exhausted batteries with new ones. Do not use batteries of different makes or not completely charged.
- This instrument is not water-tight nor resistant to splashes. Never immerse it in water or expose it to splashes
- Do not try to repair it yourself or by an unauthorized person. Always send it to Centro Style or their distributor for repair.
- The PD-meter does not generate electromagnetic disturbance which could interfere with other equipment.

This instrument is designed and manufactured for use only as explicitly stated in this manual. All other uses are considered improper.

Centro Style declines all responsibility for injury to persons or damage to property resulting from:

- Improper use and associated risks.
- Unauthorized modifications or the use of spare parts.
- uses of non original accessories.

MAINTENANCE AND CLEANING

- If not used for longer periods, remove the batteries from the instrument.
- As there are no mechanical parts, no regular maintenance is necessary. The reference lamp inside is designed for long term use.
- Clean the external surfaces with a clean cloth slightly dampened in water. Do not use solvents or potentially aggressive substances. Do not spray liquids directly onto the instrument. Do not use abrasive cleaning material.

PUTTING OUT OF SERVICE

At the end of the instrument’s working life, or if it becomes irreparably damaged, remove the batteries and dispose of them separately.

Notes of the disposal of different materials. Separation and recycling of waste material.

Refer to European Waste Disposal Directive 160200 governing the disposal of equipment or decommissioned equipment parts to dispose of or recycle.

Description of materials used in the PD-meter:

- The external structure of the body is made from ABS plastic.
- The front bearer ⑨ is made from steel and can be separated from the instrument.
- The printed circuit board with its components and wiring is inside.

TECHNICAL DATA

Power supply	4 batteries 1,5 V AA (LR 6)
Total PD range	44 – 83 mm
Partial PD range	22 – 41,5 mm
Resolution	0,5 mm
Precision	+/- 0,5 mm
Convergence distance	30 – 35 – 40 – 50 – 65 cm, 1 – 2 m, infinity
Dimensions	168 x 253 x 70 mm (including front bearer)
Weight	0,68 Kg (with batteries, without tripod)
Surrounding conditions whilst in use	Temperature between 0° e 60° C. relative humidity < 85%
CE	

WARRANTY AND SERVICE

1 CHECK THE CONTENTS OF ALL THE BOXES

2 NOTIFY CUSTOMER SERVICE OF ANY PROBLEMS WITHIN 8 DAYS

If any fault develops within the warranty period, do not open the cleaner but contact Technical Assistance immediately. All parts recognised as having manufacturing faults will be replaced under the terms of the warranty. The warranty covers faults and breakages not caused by normal wear, provided the cleaner has not been tampered with. Consumables are not covered by the warranty.

* Within one year of purchase full warranty

* from the date of invoice.

IN CASE OF MALFUNCTIONING

Never return the cleaner or parts without prior notice!
Contact your nearest authorised Service Centre before returning the cleaner or any parts, or to request repairs under the warranty.

Centro Style S.p.A. Via G.D. Martinengo, 21040 Vedano Olona, Varese.
Fax: + 39 0332 270 405
E-mail: export.dpt@centrostyle.it

Call to clarify the following points:
1 - Cleaner code and serial number
1 - Type of fault
3 - Quick over-the-phone troubleshooting and problem solving
4 - Quick over-the-phone troubleshooting with delivery of spare parts
5 - Always state the purchase date, the number of the invoice or delivery note and your customer code.

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WARNING! Read this manual carefully before using the appliance. Keep this manual safe for reference throughout the life of the appliance. Copyright Centro Style S.p.A.. All rights reserved. Total or partial reproduction of this manual is strictly forbidden without the prior written consent of Centro Style S.p.A.